# **SV38V Whole-Body Vibration Accelerometer**

The SV 38V accelerometer is dedicated for whole-body vibration measurements with the SV 106 human vibration analyser. The accelerometer has a built-in memory (TEDS) containing information about the sensitivity that is automatically transferred to the SV 106 instrument. In accordance with ISO 2631-1 seat accelerometers such as SV 38V are placed on the operator's seat either on the seat-



## **Technical Specifications**

erformance:	
Number of Axes	3
Sensitivity (± 5 %)	50 mV/(ms <sup>-2</sup> ) at 15.915 Hz
Measurement Range	0.01 ms <sup>-2</sup> RMS $\div$ 50 ms <sup>-2</sup> PEAK
Frequency Response (by design guideline, ± 3 dB)	0.1 Hz ÷ 100 Hz
Resonant Frequency	5.5 kHz (MEMS transducer)
Electrical Noise	< 0.005 ms-2 RMS, HP1 weighting
lectrical:	
iecu icai.	

Ы	lec	tri	ca	1

Supply Current	_< 5.0 mA
Supply Voltage	_5.2 V ÷ 16 V
Bias Voltage	_2.5 V ± 0.05 V
Output Impedance	_51 Ohms
Charge / Discharge Time Constant (start-up time)	_30 sec. typ.
TEDS Memory	_installed (power supply pin)

### **Environmental Conditions:**

Maximum Vibration	_100 000 ms <sup>-2</sup> shock survival for MEMS sensor
Temperature Coefficient	_<+/-0.01 %/°C
Temperature	_from -10 °C to +50 °C
Humidity	up to 90 % RH, non-condensed

#### Physical:

Sensing Element	MEMS
Cable	integrated 1.4 meters
Connector	LEMO 5-pin plug (SV 106 compatible)
Dimensions	236 mm diameter; thickness from 3.6 mm to 12 mm
Weight	550 grams (including cable and rubber cushion)

#### Accessories:

SA 38 (optional)Ca
--------------------

The policy of our company is to continually innovate and develop our products. Therefore, we reserve the right to change the specifications without prior notice.

